

ECMO Circuit Flow Impacts on General Oxygen Delivery, Synchronizing with Patient's Needs

Dr. Edgars Prozorovskis

Pauls Stradiņš Clinical University Hospital, Riga, Latvia

University of Latvia, Riga, Latvia

Riga Stradiņš University, Riga, Latvia



PAULA STRADIŅA
KLĪNISKĀ UNIVERSITĀTES
SLIMNĪCA



RĪGAS STRADIŅA
UNIVERSITĀTE



LATVIJAS
UNIVERSITĀTE

Disclosure

I have no conflict of interests to disclose

Terminology

- Veno-venous / Veno-arterial extracorporeal membrane oxygenation (**VV or VA ECMO**)
- Oxygen content (**CO₂**)
- Oxygen delivery (**DO₂**)
- Oxygen consumption (**VO₂**)
- Oxygen extraction (**VO₂**)
- Cardiac output (**CO**)
- Left ventricular outflow tract (**LVOT**)

ExtraCorporeal Life Support (ECLS)

- Mechanical devices to support heart and/or lung functions;
- A different approach to treating patients with ECLS;
- Require additional understanding of cardiopulmonary physiology, pathophysiology;
- Basic principles and physiology of the ECLS device itself

	Extracorporeal Life Support (ECLS)				
SYSTEM	Extracorporeal Membrane Oxygenation (ECMO)			Extracorporeal Carbon Dioxide Removal (ECCO ₂ R)	
SUPPORT MODE	VA ECMO	VVA ECMO	VV ECMO	VV ECCO ₂ R	AV ECCO ₂ R
CONDITION	Cardiac failure	Cardiorespiratory failure	Respiratory failure	CO ₂ retention	
APPLICATION	<ul style="list-style-type: none"> • Cardiac ECMO • ECPR • EISOR 	Cardiac and respiratory ECMO	Respiratory ECMO	Lung protection	

Patient's Needs / Cardiopulmonary Physiology

- **Normal metabolism** – amount of oxygen consumed per minute (VO_2)
 - $\text{VO}_2 = 3 \text{ ml/kg/min}$;
 - Amount of oxygen for metabolism is normally five times the amount used by the tissues ($\text{DO}_2 / \text{VO}_2 = 5:1$);
 - $\text{VO}_2 = \text{CO} \times (\text{CaO}_2 - \text{CvO}_2)$
- Oxygen content (CO_2) = **$(1.34 \times \text{Hb} \times \text{Saturation}) + (\text{PO}_2 \times 0.003)$**
 - Normally (a) $20 \text{ mlO}_2/\text{dL}$ / for C(v)O_2 $15 \text{ mlO}_2/\text{dL}$
- Oxygen delivery (DO_2) = **$\text{CO} \times \text{CaO}_2$**
 - CO (Cardiac output) measured by thermodilution or Echo: $\text{CO} = \text{HR} \times \text{SV} / \text{SV} = \text{LVOT area} \times \text{LVOT VTI}$;
 - **Normally $1000 \text{ mlO}_2/\text{min}$**


Patient's Needs / Cardiopulmonary Physiology

- DO_2 and VO_2 ratio = maintaining 5:1*

**At rest 20% of oxygen used for metabolism, 80% - leaving in the venous blood*

- Normally for patient breathing air
 - Arterial: PaO_2 – 90 mmHg / saturation – 100% / CaO_2 – 20 mLO_2/dL
 - Venous: PvO_2 – 40 mmHg / saturation – 80% / CaO_2 – 15 mLO_2/dL

Oxygen delivery (DO_2) and consumption (VO_2)



Cardiopulmonary Pathophysiology and ECMO

- **Affected DO_2 and VO_2 ratio**

- Low CO, anaemia, hypoxemia → critically low DO_2 (<2:1) → aerobic metabolism switches to anaerobic → lactate acidosis

- **ECMO**

- Cannulation strategy: size, depth/position, risk of recirculation
- Pumps: pressure
- Membrane lungs: inlet/outlet O_2 ; 100% O_2 ; size; flow rate; blood flow:gas flow = 1:1; remove of CO_2 ;
- ECMO Blood Flow was the Main Determinant of Arterial Oxygenation
- Blood flow – limitations (4-5 /min) / 60-80 ml/kg/min (dry weight) and >60% of CO
- VA ECMO – risk of Harlequin syndrome

VV ECMO – No Functional Lung

- In Severe ARDS with VV-ECMO, the Lungs Contribute Little or Nothing to Gas Exchange
- Right atrium (RA): oxygenated blood from ECMO + venous blood from systemic organs
 - Arterial saturation ranges from 60 to 90 % (*depending on ECMO flow, native venous flow, lung function and cardiac output*)
 - **NORMAL SYSTEMATIC OXYGEN DELIVERY** (if CO and Hbg conc. adequate)
 - Hypoxemia (PaO₂ 40-60 mmHg, SaO₂ 70-90%) always occurs with VV support and is adequate to maintain normal oxygen delivery
- CO = ECMO flow + Native venous flow

Date for DO_2/VO_2 Consumption Ratio

- **Patient monitoring**
 - Ventilation parameters
 - Hemodynamic
 - Arterial and venous blood gases analyse
 - Haemoglobin/ haematocrit levels
 - Cardiac function
- **ECMO system monitoring**
 - Flow
 - Pressures
 - Inlet-outlet oxygen content

Calculation of DO₂/VO₂ Ratio During ECMO

Step 1 Calculate Oxygen Content	Arterial Blood Oxygen Content <small>(mL O₂/L)</small>	=	1.34 <small>(mL O₂/g Hb)</small>	x	Hemoglobin <small>(g/dL)</small>	x	$\frac{10 \text{ dL}}{1 \text{ L}}$ <small>unit conversion</small>	x	Arterial SaO ₂ <small>(decimal)</small>	+	$\left(\text{Arterial pO}_2 \times 0.0031 \times \frac{10 \text{ dL}}{1 \text{ L}} \right)$ <small>(mm Hg) (mL O₂/dL/mm Hg) unit conversion</small>
	Pre-Oxygenator Blood Oxygen Content <small>(mL O₂/L)</small>	=	1.34 <small>(mL O₂/g Hb)</small>	x	Hemoglobin <small>(g/dL)</small>	x	$\frac{10 \text{ dL}}{1 \text{ L}}$ <small>unit conversion</small>	x	Pre-Oxygenator SaO ₂ <small>(decimal)</small>	+	$\left(\text{Pre-Oxygenator pO}_2 \times 0.0031 \times \frac{10 \text{ dL}}{1 \text{ L}} \right)$ <small>(mm Hg) (mL O₂/dL/mm Hg) unit conversion</small>
	Post-Oxygenator Blood Oxygen Content <small>(mL O₂/L)</small>	=	1.34 <small>(mL O₂/g Hb)</small>	x	Hemoglobin <small>(g/dL)</small>	x	$\frac{10 \text{ dL}}{1 \text{ L}}$ <small>unit conversion</small>	x	Post-Oxygenator SaO ₂ <small>(decimal)</small>	+	$\left(\text{Post-Oxygenator pO}_2 \times 0.0031 \times \frac{10 \text{ dL}}{1 \text{ L}} \right)$ <small>(mm Hg) (mL O₂/dL/mm Hg) unit conversion</small>

- Constant 13.4 mL O₂/g Hb: accounts for the fact that 1.34 ml of O₂ is carried per g of Hb (13.4 is used in the equation to correct the units from dL to L)
- Constant 0.0031 mL O₂/L/mm Hg: solubility coefficient of oxygen at body temperature

Calculation of DO₂/VO₂ Ratio During ECMO

Step 2 Calculate Total Flow	$\text{Total Flow} = \text{VV-ECMO Flow} \times \frac{\text{Post-Oxygenator Blood Oxygen Content} - \text{Pre-Oxygenator Blood Oxygen Content}}{\text{Arterial Blood Oxygen Content} - \text{Pre-Oxygenator Blood Oxygen Content}}$	
	<table border="1" style="width: 100%;"> <tr> <td style="text-align: center; width: 50%;"> VV-ECMO Flow Rate as Set on the Machine (L/min) </td> <td style="text-align: center; width: 50%;"> $\text{Venous Flow (L/min)} = \text{Total Flow (L/min)} - \text{VV-ECMO Flow (L/min)}$ </td> </tr> </table>	VV-ECMO Flow Rate as Set on the Machine (L/min)
VV-ECMO Flow Rate as Set on the Machine (L/min)	$\text{Venous Flow (L/min)} = \text{Total Flow (L/min)} - \text{VV-ECMO Flow (L/min)}$	

- When Two Blood Flows Containing Different Oxygen Contents Mix, the Resultant Oxygen Content is the Average of the Amount of Oxygen Content in Each of the Two Flows
- **Total oxygen compsuption:** $[(\text{Post-Oxygenator Oxygen Content} \times \text{VV-ECMO Flow}) / \text{Total Flow}] + [(\text{Pre-Oxygenator Oxygen Content} \times \text{Native Venous Flow}) / \text{Total Flow}]$

Calculation of DO₂/VO₂ Ratio During ECMO

Step 3
Calculate Oxygen Delivery/Consumption Ratio

Oxygen Delivery (mL O ₂ /min)	=	Arterial Blood Oxygen Content (mL O ₂ /L)	x	Total Flow (L/min)
Oxygen Consumption (mL O ₂ /min)	=	Weight (kg)	x	3 (mL O ₂ /min/kg)

Oxygen Delivery/Consumption Ratio = $\frac{\text{Oxygen Delivery}}{\text{Oxygen Consumption}}$

- Oxygen Consumption = Weight (in kg) x 3 mL O₂/kg/min

Maintain Oxygen Delivery/Oxygen Consumption Ratio >3

Sample Calculation of Oxygen Delivery (DO₂)/Oxygen Consumption (VO₂) Ratio During Venovenous Extracorporeal Membrane Oxygenation (VV-ECMO)

Step 1
Calculate Oxygen Content

$$\text{Arterial Blood Oxygen Content (mL O}_2\text{/L)} = 1.34 \times 10.5 \times \frac{10 \text{ dL}}{1 \text{ L}} \times 0.88 + \left(50 \times 0.0031 \times \frac{10 \text{ dL}}{1 \text{ L}} \right) = 125.4$$

(mL O₂/L) (mL O₂/g Hb) (g/dL) unit conversion (decimal) (mm Hg) (mL O₂/dL/mm Hg) unit conversion (mL O₂/L)

$$\text{Pre-Oxygenator Blood Oxygen Content (mL O}_2\text{/L)} = 1.34 \times 10.5 \times \frac{10 \text{ dL}}{1 \text{ L}} \times 0.64 + \left(40 \times 0.0031 \times \frac{10 \text{ dL}}{1 \text{ L}} \right) = 91.3$$

(mL O₂/L) (mL O₂/g Hb) (g/dL) unit conversion (decimal) (mm Hg) (mL O₂/dL/mm Hg) unit conversion (mL O₂/L)

$$\text{Post-Oxygenator Blood Oxygen Content (mL O}_2\text{/L)} = 1.34 \times 10.5 \times \frac{10 \text{ dL}}{1 \text{ L}} \times 1.00 + \left(300 \times 0.0031 \times \frac{10 \text{ dL}}{1 \text{ L}} \right) = 150.0$$

(mL O₂/L) (mL O₂/g Hb) (g/dL) unit conversion (decimal) (mm Hg) (mL O₂/dL/mm Hg) unit conversion (mL O₂/L)

Step 2
Calculate Total Flow

$$\text{Total Flow} = 4 \text{ L/min} \times \frac{150.0 - 91.3}{125.4 - 91.3} = 6.9 \text{ (L/min)}$$

VV-ECMO Flow Rate of 4L/min
(L/min)

$$\text{Venous Flow} = \text{Total Flow} - \text{VV-ECMO Flow} = 2.9 \text{ (L/min)}$$

Step 3
Calculate Oxygen Delivery/Consumption Ratio

$$\text{Oxygen Delivery (mL O}_2\text{/min)} = 125.4 \times 6.9 = 865.3$$


(mL O₂/min) (mL O₂/L) (L/min)

$$\text{Oxygen Consumption (mL O}_2\text{/min)} = 70 \times 3 = 210.0$$

(mL O₂/min) (kg) (mL O₂/min/kg)

$$\text{Oxygen Delivery/Consumption Ratio} = \frac{865.3}{210.0} = 4.1$$

Mixing of Venovenous Extracorporeal Membrane Oxygenation (VV-ECMO) Blood Flow with Native Venous Blood Flow



Example from Our Department

Case Nr.	N18	
ECMO day	5	
Patient date		
Age (y)	20	
Weight (kg)	70	
Height (cm)	160	
ECMO support	VV	
Cannulas		
Out.	v.femoralis dx	26 Fr
In.	v.jug.int.dx	19 Fr
Distance (cm)	10	

Oxygen content		
CaO ₂ (arterial)	133	mLO ₂ /L
CaO ₂ (pre)	113	mLO ₂ /L
CaO ₂ (post)	146	mLO ₂ /L

ABG					
Arterial blood		Pre-oxygenator blood		Post-oxygenator blood	
pO ₂ (mmHg)	91	pO ₂ (mmHg)	50	pO ₂ (mmHg)	469
SO ₂ (decimal)	0,98	SO ₂ (decimal)	0,84	SO ₂ (decimal)	0,99

Hemoglobin (g/dL)	9,9
-------------------	-----

Cardiac output (Echo) =	0
-------------------------	---

HR (x')	
LVOT AREA	
LVOT VTI	


Constants	
1,34	mLO ₂ /gHb
0,0031	mLO ₂ /L/mmHg
3	mLO ₂ /kg/min
1	L
10	dL

Flows		
ECMO flow	4,0	l/min
Venous flow	2,6	l/min
Total flow	6,6	l/min

Oxygen delivery / consumption ratio	
DO ₂ (mLO ₂ /min)	603
VO ₂ (mLO ₂ /min)	210
DO₂/VO₂ ratio	2,9

ECMO Management for Patients Needs

- **Ventilation setup**
 - NO IDEAL setups, only recommendations and protective strategies
 - Turning up the ventilator FiO₂ or airway pressure will not help
- **↑ Blood oxygen content**
 - ↑ ECMO flow rate
 - ↑ haemoglobin/ haematocrit (up to 12 – 14 g/dL)
- **↓ Recirculation**
 - ↑ distance between cannulas
- **↓ Oxygen consumption**
 - ↑ sedation
 - Therapeutic hypothermia
 - Neuromuscular blockage
- **Cardiac output / intrapulmonary shunt**
 - B-blockers
- *Switch to venoarterial or hybrid configuration*



Take Home Messages

- **VO₂/DO₂ should be calculated for every patients on ECMO**
- Adjust the ECMO variables and the patient variables to maintain DO₂/VO₂ =/> 3:1;
- Non damaging ventilation: (FiO₂ 0.3 – 0.5, CPAP 15-20cm H₂O, RR 5 – 10, Pplat<25);
- Titrated ECMO flow to maintain SaO₂ 90%;
- Adjust the sweep gas to keep the PaCO₂ 40 mmHg;
- If low SaO₂ with adequate ECMO flow -> recirculation/oxygenator dysfunction or high VO₂ (↑ ECMO flow? ↑Hb? Temp.? ↓ CO ? – which one?)
- **Allow healing of the lungs and rest**

Take Home Messages

- When the patient is stable (*usually 6-12 hours*) determine the variables of O₂ kinetics:
 - If oxygen supply is adequate (DO₂:VO₂ over 3) no changes are necessary;
 - If oxygen supply is inadequate (DO₂:VO₂ under 3) and the patient is anaemic, transfuse to a higher haemoglobin;
 - If DO₂ is still inadequate change the drainage cannula to a larger size and increase flow
- Remember:
 - DO₂ ↑, if ↑ Hgb / CO / PaO₂
 - VO₂ ↓, if – hypothermia, anaesthesia, decrease work of breathing

Thank You!



References

1. Rosen, I.M. and Manaker S. (2023). Oxygen delivery and consumption. *UpToDate*. Retrieved September 24, 2023 from <https://www.uptodate.com/contents/oxygen-delivery-and-consumption>
2. Venovenous Extracorporeal Membrane Oxygenation (VV-ECMO). MD Nexus. Retrieved September 24, 2023 from <https://www.mdnx.com/topics-2/procedures/venovenous-extracorporeal-membrane-oxygenation/>
3. Conrad, S.A. *et al.* (2018) “The Extracorporeal Life Support Organization Maastricht Treaty for Nomenclature in Extracorporeal Life Support A Position Paper of the Extracorporeal Life Support Organization,” *American journal of respiratory and critical care medicine*, 198(4), pp. 447–451.
4. Bartlett, R.H. (2023). Physiology of Extracorporeal Gas Exchange. In *Comprehensive Physiology*, Y.S. Prakash (Ed.). <https://doi.org/10.1002/cphy.c190006>
5. Tonna, Joseph E *et al.* “Management of Adult Patients Supported with Venovenous Extracorporeal Membrane Oxygenation (VV ECMO): Guideline from the Extracorporeal Life Support Organization (ELSO).” *ASAIO journal (American Society for Artificial Internal Organs : 1992)* vol. 67,6 (2021): 601-610.
6. Marhong, J.D. *et al.* (2014) “Mechanical ventilation during extracorporeal membrane oxygenation an international survey,” *Annals of the American Thoracic Society*, 11(6), pp. 939–944.
7. Togo, K. *et al.* (2015) “Impact of bypass flow rate and catheter position in veno-venous extracorporeal membrane oxygenation on gas exchange in vivo,” *Journal of artificial organs*, 18(2), pp. 128–135.